*Introduction*

*Einführung*

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*Introduction*

*Introdução*

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## B52 / B53

Mounting Instruction

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
## **Life-Size Dual-Sex 3B Scientific® Muscle Figure B52 Asian / B53 European**

This high-quality 3B Scientific® Muscle Figure was developed especially for training and studies at university level to meet the high demands of a good study of anatomy.

The 3B Scientific® Muscle Figure includes:

### **CD ROM - 3B Scientific® Muscle Figure**

19 color illustrations of all structures and organs with the respective explanations in 7 different languages (Latin, English, German, Spanish, Portuguese, French, Japanese) in pdf file format.

### **CD ROM - 3B Scientific® ANATOMYtrainer™**

Interactive learning program with quiz mode for students of medicine. With 400 digital illustrations and 3,000 relevant terms needed to prepare for examinations. In 6 different languages (Latin, English, German, Spanish, Portuguese, French).

### **CD ROM - 3B Scientific® MUSCLEtrainer™**

Interactive learning program with quiz mode for students of medicine. With 248 digital illustrations and over 400 relevant terms needed to prepare for examinations. In 6 different languages (Latin, English, German, Spanish, Portuguese, French).

### **How to assemble the 3B Scientific® Muscle Figure**

To ensure that even smallest anatomical details are represented as accurately as possible, we chose to manufacture this figure out of a special high-quality plastic. This material is a little more sensitive compared to ordinary materials. We therefore ask you to proceed carefully during the assembly and disassembly of the figure and follow these steps:

Work with a partner when assembling the figure.

Unpack the following partially assembled model parts:

- Torso body with skin arm and leg (includes internal organs, chest wall and neck muscle)
- Muscular upper leg with gluteal muscle
- Muscular lower leg
- Muscle arm
- Head
- 1 genital insert
- Base with washer and wing nut

- Cover a table with a soft cloth and place the torso body with skin arm and skin leg onto the table. The skin leg should protrude slightly (approx. 10-20 cm) from the table edge.
- Slide the base with the hole over the threaded bolt protruding from the foot and fasten it with the washer and wing nut.
- 2 of the 4 rollers under the base are equipped with brakes. Please lock these during assembly.
- Carefully raise the figure.
- Connect the muscular upper leg (without gluteal muscle) with the muscular lower leg and fasten both retaining hooks at the knee. Attach the assembled muscular leg to the pelvis using the guide pin.
- Now insert the gluteal muscle into the guide hole from the top.
- Attach the muscular arm to the shoulder using the guide pin and fasten the retaining hook.
- Place the head onto the neck using both connecting pins.
- Remove the wrapping material inserted between the individual model parts in order to protect the painted surfaces during transport.

### **Handling and Care of the 3B Scientific® Muscle Figure**

- The arm and leg muscles are equipped with self-locking connectors and can be removed in the correct anatomical order.
- To remove the muscular leg it is always necessary to first remove the gluteal muscle.
- To remove the chest wall it is always necessary to first remove the neck muscle. The chest wall should be tilted forwards slightly and then pulled off in horizontal direction.
- The internal organs can be removed in their correct anatomical order.
- To exchange the genital inserts please remove the internal organs first.
- Please ensure that the figure is completely reassembled after each demonstration.
- Avoid strong thermal stress such as prolonged direct exposure to sun and only use mild soap solutions to clean the figure (no solvents!)
- All organs and assembly parts can be reordered upon exact specification. Thanks to 3B Scientific high standard of quality, however, we are sure that you will only make use of this opportunity when parts are lost.

And now enjoy your 3B Scientific® Muscle Figure.

Further information on 3B Scientific® products and many new items is available on the Internet at:  
<http://www.3bscientific.com>

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## Introduction

This pdf file is intended to assist you in preparing and holding your lectures. It identifies all structures and organs of the 3B Scientific® Muscle Figure B52/B53 on a total of 19 color illustrations and in 7 different languages. In addition, you will find many informative background texts for your lessons.

Of course you can also make color printouts of specific pages for your students or print your own overhead sheets.

### And this is how it works:

- Simply click on one of the blue fields with page numbers and titles (e.g. #15 Heart) on the right hand side of this page and you will automatically go to the relevant color illustration. The numbered reference lines pointing to the individual anatomic structures allow for quick orientation in the text section.
- Left of the illustrations you will find national flags representing the languages Latin, English, German, Spanish, French, Portuguese and Japanese. Click on your language of choice to display the text of the structures numbered in the illustration in this language.
- To return to the previous page simply click on the ◀ symbol in the menu bar.
- Opening the [navigation window](#) additionally allows you to navigate through a directory tree.

- In case you have lost the supplied mounting instructions, simply click here:

[Mounting Instruction B52/B53](#)

- #01 Trunk bisexual
- #02 Trunk, with open back section
- #03 Hip and thigh
- #04 Muscular arm
- #05 Muscular arm
- #06 Muscular leg
- #07 Muscular leg
- #08 Musclar lower leg with knee
- #09 Musclar lower leg with knee
- #10 Chest wall and abdominal wall muscles
- #11 Head
- #12 Head with neck
- #13 Brain hemisphere
- #14 Lungs
- #15 Heart
- #16 Liver and open stomach
- #17 Small intestine and large intestine
- #18 Male genital organs
- #19 Female genital organs

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



## Muscular arm

- |    |  |    |   |
|----|--|----|---|
| 1  | Tendinous junction                           | 48 | Lumbrical muscles of hand                       |
| 2  | Dorsal interosseous muscles of hand          | 49 | Tendon of flexor digitorum superficialis muscle |
| 3  | Radial nerve                                 | 50 | Superficial transverse metacarpal ligament      |
| 4  | Radial artery                                | 51 | Common palmar digital arteries                  |
| 5  | Extensor retinaculum                         | 52 | Digital synovial sheath                         |
| 6  | Abductor pollicis brevis muscle              | 53 | Tendon of flexor digitorum profundus muscle     |
| 7  | Extensor carpi radialis brevis muscle        | 54 | Annular part of fibrous tendon sheaths          |
| 8  | Extensor carpi radialis longus muscle        | 55 | Proper palmar digital nerves of ulnar nerve     |
| 9  | Brachioradialis muscle                       | 56 | Cruciform part of fibrous tendon sheaths        |
| 10 | Brachialis muscle                            | 57 | Superficial palmar arch                         |
| 11 | Deltoid muscle                               | 58 | Abductor digiti minimi muscle of hand           |
| 12 | Supraspinatus muscle                         | 59 | Opponens digiti minimi muscle of hand           |
| 13 | Scapular spine                               | 60 | Superficial branch of ulnar nerve               |
| 14 | Trapezius muscle                             | 61 | Dorsal metatarsal arteries                      |
| 15 | Infraspinatus muscle                         | 62 | Head of metacarpal bone                         |
| 16 | Latissimus dorsi muscle                      | 63 | Dorsal digital arteries of hand                 |
| 17 | Teres major muscle                           | 64 | Dorsal digital nerves of fingers                |
| 18 | Teres minor muscle                           | 65 | Tendon of extensor indicis muscle               |
| 19 | Triangular space                             | 66 | Body of metacarpal bone                         |
| 20 | Long head of triceps brachii muscle          | 67 | Tendon of extensor pollicis muscle              |
| 21 | Lateral head of triceps brachii muscle       |    |   |
| 22 | Tendon of triceps brachii muscle             |    |   |
| 23 | Ulnar nerve                                  |    |   |
| 24 | Anconeus muscle                              |    |   |
| 25 | Medial epicondyle of humerus                 |    |   |
| 26 | Olecranon process of ulna                    |    |   |
| 27 | Flexor carpi ulnaris muscle                  |    |   |
| 28 | Extensor digitorum muscle of hand            |    |   |
| 29 | Extensor carpi ulnaris muscle                |    |   |
| 30 | Tendon of extensor digitorum muscle of hand  |    |   |
| 31 | Head of ulna                                 |    |   |
| 32 | Tendon of extensor digiti minimi muscle      |    |   |
| 33 | Dorsal branch of ulnar nerve                 |    |   |
| 34 | Coracobrachialis muscle                      |    |   |
| 35 | Axillary nerve                               |    |   |
| 36 | Posterior circumflex humeral artery          |    |   |
| 37 | Quadrangular space                           |    |   |
| 38 | Brachial artery                              |    |   |
| 39 | Humerus                                      |    |   |
| 40 | Lateral epicondyle of humerus                |    |   |
| 41 | Tendon of flexor carpi radialis muscle       |    |   |
| 42 | Median nerve                                 |    |   |
| 43 | Flexor retinaculum                           |    |   |
| 44 | Tendon of palmaris longus muscle             |    |   |
| 45 | Common palmar digital nerves of median nerve |    |   |
| 46 | Flexor pollicis brevis muscle                |    |   |
| 47 | Adductor pollicis muscle                     |    |   |

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


## Muscular arm

- 1 Trapezius muscle
- 2 Clavicle
- 3 Subclavius muscle
- 4 Axillary artery
- 5 Pectoralis minor muscle
- 6 Brachial plexus
- 7 Deltoideus muscle
- 8 Pectoralis major muscle
- 9 Brachial artery
- 10 Median nerve
- 11 Ulnar nerve
- 12 Short [medial] head of biceps brachii muscle
- 13 Tendon of biceps brachii muscle
- 14 Brachioradialis muscle
- 15 Flexor digitorum superficialis muscle
- 16 Ulnar artery
- 17 Flexor retinaculum
- 18 Abductor digiti minimi muscle of hand
- 19 Head of ulna
- 20 Extensor carpi ulnaris muscle
- 21 Flexor carpi ulnaris muscle
- 22 Flexor carpi radialis muscle
- 23 Pronator teres muscle
- 24 Palmaris longus muscle
- 25 Medial epicondyle of humerus
- 26 Brachialis muscle
- 27 Medial head of triceps brachii muscle
- 28 Long head of triceps brachii muscle
- 29 Latissimus dorsi muscle
- 30 Teres major muscle
- 31 Thoracodorsal nerve
- 32 Subscapularis muscle
- 33 Medial margin of scapula
- 34 Coracobrachialis muscle
- 35 Flexor digitorum profundus muscle
- 36 Dorsal branch of ulnar nerve
- 37 Ulna
- 38 Olecranon process of ulna
- 39 Inferior collateral ulnar artery
- 40 Humerus
- 41 Triceps brachii muscle
- 42 Radial nerve
- 43 Acromion [acromial bone]
- 44 Musculocutaneous nerve
- 45 Radial artery
- 46 Anterior interosseous artery
- 47 Pronator quadratus muscle
- 48 Tendon of flexor carpi radialis muscle
- 49 Abductor pollicis brevis muscle
- 50 Superficial palmar arch
- 51 Tendon of palmaris longus muscle
- 52 Tendon of flexor pollicis longus muscle

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



## Muscular leg

- 1 Trapezius muscle
- 2 Clavicle
- 3 Subclavius muscle
- 4 Axillary artery
- 5 Pectoralis minor muscle
- 6 Brachial plexus
- 7 Deltoideus muscle
- 8 Pectoralis major muscle
- 9 Brachial artery
- 10 Median nerve
- 11 Ulnar nerve
- 12 Short [medial] head of biceps brachii muscle
- 13 Tendon of biceps brachii muscle
- 14 Brachioradialis muscle
- 15 Flexor digitorum superficialis muscle
- 16 Ulnar artery
- 17 Flexor retinaculum
- 18 Abductor digiti minimi muscle of hand
- 19 Head of ulna
- 20 Extensor carpi ulnaris muscle
- 21 Flexor carpi ulnaris muscle
- 22 Flexor carpi radialis muscle
- 23 Pronator teres muscle
- 24 Palmaris longus muscle
- 25 Medial epicondyle of humerus
- 26 Brachialis muscle
- 27 Medial head of triceps brachii muscle
- 28 Long head of triceps brachii muscle
- 29 Latissimus dorsi muscle
- 30 Teres major muscle
- 31 Thoracodorsal nerve
- 32 Subscapularis muscle
- 33 Medial margin of scapula
- 34 Coracobrachialis muscle
- 35 Flexor digitorum profundus muscle
- 36 Dorsal branch of ulnar nerve
- 37 Ulna
- 38 Olecranon process of ulna
- 39 Inferior collateral ulnar artery
- 40 Humerus
- 41 Triceps brachii muscle
- 42 Radial nerve
- 43 Acromion [acromial bone]
- 44 Musculocutaneous nerve
- 45 Radial artery
- 46 Anterior interosseous artery
- 47 Pronator quadratus muscle
- 48 Tendon of flexor carpi radialis muscle
- 49 Abductor pollicis brevis muscle
- 50 Superficial palmar arch
- 51 Tendon of palmaris longus muscle
- 52 Tendon of flexor pollicis longus muscle
- 53 Quadratus femoris muscle

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
## Musclar lower leg with knee

- |   |  |
|---|--|
| 1 Patellar tendon                               | 48 Peroneus longus muscle                            |
| 2 Adipose body of patella                       | 49 Extensor digitorum longus muscle (of foot)        |
| 3 Medial patellar retinaculum                   | 50 Lateral malleolus                                 |
| 4 Transverse ligament of knee                   | 51 Transverse ligament of ankle                      |
| 5 Anterior cruciate ligament                    | 52 Flexor digitorum longus muscle                    |
| 6 Synovial membrane                             | 53 Soleus muscle                                     |
| 7 Medial meniscus of knee                       | 54 Tibia [shin bone]                                 |
| 8 Sartorius muscle                              | 55 Tuberosity of tibia                               |
| 9 Medial condyle of tibia                       | 56 Tendon of sartorius muscle                        |
| 10 Tendon of gracilis muscle                    | 57 Extensor digitorum brevis muscle of foot          |
| 11 Posterior cruciate ligament                  | 58 Extensor hallucis brevis muscle                   |
| 12 Tendon of semitendinosus muscle              | 59 Cruciate ligament of ankle                        |
| 13 Tendon of semimembranosus muscle             | 60 Tendon of peroneus brevis muscle                  |
| 14 Popliteal artery                             | 61 Tendon of peroneus tertius muscle                 |
| 15 Popliteal vein                               | 62 Dorsal interosseous muscles of foot               |
| 16 Tibial nerve                                 | 63 Tendons of extensor digitorum longus pedis muscle |
| 17 Semimembranosus muscle                       | 64 Superficial transverse metatarsal ligament        |
| 18 Lateral condyle of tibia                     | 65 Cruciform part of fibrous tendon sheaths          |
| 19 Lateral meniscus of knee                     | 66 Annular part of fibrous tendon sheaths            |
| 20 Long head of biceps femoris muscle           | 67 Superficial transverse metacarpal ligament        |
| 21 Tendon of popliteal muscle                   | 68 Tendon of flexor digitorum brevis muscle          |
| 22 Short head of biceps femoris muscle          | 69 Tendon of flexor digitorum longus muscle          |
| 23 Subpopliteal recess                          | 70 Flexor digiti minimi brevis muscle of foot        |
| 24 Lateral patellar retinaculum                 | 71 Abductor digiti minimi muscle of foot             |
| 25 Medial head of gastrocnemius muscle          | 72 Flexor digitorum brevis muscle                    |
| 26 Medial condyle of femur                      | 73 Calcaneal tuber                                   |
| 27 Femur [thigh bone]                           | 74 Plantar aponeurosis                               |
| 28 Patellar surface of femur                    | 75 Abductor hallucis muscle                          |
| 29 Deep infrapatellar bursa                     | 76 Tendon of flexor hallucis longus muscle           |
| 30 Infrapatellar synovial fold                  | 77 Flexor hallucis brevis muscle                     |
| 31 Patellar ligament                            | 78 Lumbrical muscle of foot                          |
| 32 Common fibular nerve                         | 79 Plantaris muscle                                  |
| 33 Lateral condyle of femur                     |  |
| 34 Biceps femoris muscle                        |  |
| 35 Iliotibial tract                             |  |
| 36 Vastus lateralis muscle                      |  |
| 37 Tendon of rectus femoris muscle              |  |
| 38 Vastus intermedius muscle                    |  |
| 39 Vastus medialis muscle                       |  |
| 40 Tendon of adductor magnus muscle             |  |
| 41 Saphenous nerve                              |  |
| 42 Gracilis muscle                              |  |
| 43 Semitendinosus muscle                        |  |
| 44 Patella                                      |  |
| 45 Head of fibula                               |  |
| 46 Lateral [great] head of gastrocnemius muscle |  |
| 47 Tibialis anterior muscle                     |  |

\* This figure shows a section that is only visible in the partial model M22 (muscular lower leg with knee).

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


## Head with neck

- |  |  |
|--|--|
| 1 Medial branch of supraorbital nerve            | 48 Zygomatic and buccal branches of facial nerve |
| 2 Lateral branch of supraorbital nerve           | 49 Lateral crus of major alar cartilage          |
| 3 Frontal belly of occipitofrontalis muscle      | 50 Temporal branches of facial nerve             |
| 4 Galea aponeurotica                             | 51 External carotid artery                       |
| 5 Temporalis muscle                              | 52 Lingual tonsil                                |
| 6 Palpebral part of orbicularis oculi muscle     | 53 Palatine tonsil                               |
| 7 Orbital part of orbicularis oculi muscle       | 54 Hyoglossus muscle                             |
| 8 Superficial temporal artery                    | 55 Constrictor pharyngis superior muscle         |
| 9 Nasal muscle                                   | 56 Greater horn of hyoid bone                    |
| 10 Levator labii superioris alaeque nasi muscle  | 57 Internal carotid artery                       |
| 11 Levator labii superioris muscle               | 58 Scalenus medius muscle                        |
| 12 External ear                                  | 59 Cervical nerves II-IV                         |
| 13 Zygomaticus minor muscle                      | 60 Accessory nerve                               |
| 14 Zygomaticus major muscle                      | 61 Brachial plexus                               |
| 15 Parotid duct [duct of STENON]                 | 62 Scalenus anterior muscle                      |
| 16 Parotid gland                                 | 63 Internal jugular vein                         |
| 17 Zygomatic arch                                | 64 Common carotid artery                         |
| 18 Masseter muscle                               | 65 Constrictor pharyngis inferior muscle         |
| 19 Buccinator muscle                             | 66 Thyroid membrane                              |
| 20 Risorius muscle                               | 67 Thyrohyoideus muscle                          |
| 21 Facial artery                                 | 68 Stylohyoideus muscle                          |
| 22 Mandibula [mandible]                          | 69 Thyroid cartilage                             |
| 23 Depressor anguli oris muscle                  | 70 Laryngeal prominence                          |
| 24 Inferior labial artery                        | 71 Sternohyoideus muscle                         |
| 25 Depressor labii inferioris muscle             | 72 Medial thyrohyoid ligament                    |
| 26 Superior labial artery                        | 73 Body of hyoid bone                            |
| 27 Mental muscle                                 | 74 Mylohyoideus muscle                           |
| 28 Orbicularis oris muscle                       | 75 Genohyoideus muscle                           |
| 29 Greater alar cartilage                        | 76 Genioglossus muscle                           |
| 30 Angular artery                                |  |
| 31 Frontal branch of superficial temporal artery |  |
| 32 Transverse facial artery                      |  |
| 33 Zygomaticoorbital artery                      |  |
| 34 Temporoparietalis muscle                      |  |
| 35 Anterior auricular artery                     |  |
| 36 Auricularis superior muscle                   |  |
| 37 Occipital artery                              |  |
| 38 Auriculotemporal nerve                        |  |
| 39 Greater occipital nerve                       |  |
| 40 Auricularis anterior muscle                   |  |
| 41 Auricularis posterior muscle                  |  |
| 42 Posterior auricular artery                    |  |
| 43 Lesser occipital nerve                        |  |
| 44 Trapezius muscle                              |  |
| 45 Levator scapulae muscle                       |  |
| 46 Sternocleidomastoideus muscle                 |  |
| 47 Great auricular nerve                         |  |

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



The heart is the central circulatory organ and functions like a pump. It is a hollow muscular body which contracts 60 times (at rest) to 180 times (under stress) per minute. Approx. 70ml of blood is ejected with every one of these contractions. A septum divides the heart into a left and a right section. Each of these areas has an atrium and a ventricle. In the right section of the heart, deoxygenated blood from the systemic circulation is collected and pumped into the lungs to be reoxygenated. The blood from the pulmonary circulation then enters the left section of the heart and from there is then pumped to all parts of the body. A blood corpuscle takes approx. one minute to cover these distances. Valves to prevent reflux of blood may be found between the different sections of the heart.

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



## Liver and open stomach

### a Liver

- 1 Bare area of liver
- 2 Diaphragmatic surface of hepar
- 3 Visceral surface of hepar
- 4 Colic impression of liver
- 5 Duodenal impression of liver
- 6 Gastric impression of liver
- 7 Esophageal impression of liver
- 8 Renal impression of liver
- 9 Suprarenal impression of liver
- 10 Caudate lobe of liver
- 11 Right lobe of liver
- 12 Quadrate lobe of liver
- 13 Left lobe of liver
- 14 Inferior margine of liver
- 15 Caudate process of liver
- 16 Common hepatic artery
- 17 Right hepatic duct
- 18 Falciform ligament of liver
- 19 Left hepatic duct
- 20 Inferior vena cava
- 21 Hepatic veins
- 22 Neck of gallbladder
- 23 Body of gallbladder
- 24 Fundus of gallbladder
- 25 Cystic duct
- 26 Common bile duct
- 27 Portal vein of liver

### b Stomach

- 1 Right gastric artery
- 2 Left gastric artery
- 3 Short gastric arteries
- 4 Right gastroepiploic artery
- 5 Left gastroepiploic artery
- 6 Major gastric curvature
- 7 Minor gastric curvature
- 8 Angular notch of stomach
- 9 Cardiac notch of stomach
- 10 Cardiac part of stomach
- 11 Body of stomach
- 12 Fundus of stomach
- 13 Abdominal part of esophagus
- 14 Pylorus
- 15 Mucosa of stomach
- 16 Muscular coat of stomach
- 17 Gastric folds
- 18 Pyloric antrum [pyloric canal]
- 19 Greater omentum

### The liver

The liver is the largest organ in the human body and weighs approx. 1.5kg in the adult. It may be termed a "central laboratory" in terms of its tasks in the organism. Numerous metabolic processes take place here, including

- protein synthesis and breakdown
- synthesis, breakdown and storage of glycogen, the energy carrier in carbohydrate metabolism
- synthesis and breakdown of fatty acids
- detoxification of foreign substances such as alcohol, drugs and medicines
- hormone inactivation
- production of approx. 1L bile per day. This is thickened and stored in the gallbladder. Its tasks include lipolysis and excretion of toxic metabolic products.

The upper side of the liver is attached to the diaphragm, its base is in contact with the stomach and intestine.

### The stomach

Food taken up by the mouth and oesophagus is stored temporarily in the stomach and passed on in small amounts to the intestine. The cells of the stomach wall produce hydrochloric acid to kill bacteria taken in with the food and various enzymes start to partly digest the food. The partly digested food passes through the pylorus and enters the duodenum.

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



## Small intestine and large intestine

- 1 Pancreatic duct
- 2 Ascending colon
- 3 Descending colon
- 4 Sigmoid colon [pelvic colon]
- 5 Transverse colon
- 6 Haustra of colon [Sacculations of colon]
- 7 Ostium of ileocecal valve
- 8 Ileocecal valve
- 9 Pancreatic head
- 10 Pancreatic cauda [pankreatic tail]
- 11 Pancreatic body
- 12 Pancreas
- 13 Duodenum
- 14 Inferior duodenal flexure
- 15 Superior duodenal flexure
- 16 Ascending part of duodenum
- 17 Descending part of duodenum
- 18 Common bile duct
- 19 Splenic artery and vein [lienial artery and vein]
- 20 Vermiform appendage
- 21 Ostium of vermiform appendage
- 22 Cecum [blind sac]
- 23 Free taenia
- 24 Ileum
- 25 Jejunum
- 26 SANTORINI's duct
- 27 Superior mesenteric artery and vein
- 28 Inferior mesenteric artery and vein
- 29 Jejunal and ileal arteries and veins
- 30 Right colic artery and vein
- 31 Middle colic artery and vein
- 32 Left colic artery and vein
- 33 Sigmoid arteries and veins
- 34 Superior rectal artery
- 35 Mesentery of transverse colon

The intestine is made up of the small and large intestine and the first part of the small intestine is called the duodenum. The excretory duct of the gallbladder merges with the pancreatic duct which opens into the duodenum at the major duodenal papilla. The further parts of the small intestine are the jejunum and the ileum. The total length of the small intestine in living humans cannot be stated exactly due to constant contractions and movements. In a relaxed state, the small intestine is approx. 5 - 6m long. If it is cut open, numerous mucosal folds may be seen inside which also feature villi which give the mucosa a velvety appearance. The folds and villi increase the surface of the small intestine to 5 - 10 m<sup>2</sup>. This large surface is used for effective absorption of nutrients. The digestive juices from the liver and pancreas added in the duodenum first split up food into molecules which can be absorbed into the blood and lymph vessels in the further course of the small intestine. Muscles in the intestinal wall produce contraction waves which propel the chyme forwards. Finally, the remaining chyme is thickened by water absorption and excreted in the large intestine which is approx. 1.5m long.

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



## Male genital organs

- |   |                    |
|---|--------------------|
| 1 Rectus abdominis muscle                               | 48 Pubic symphysis |
| 2 Sigmoid colon [pelvic colon]                          |                    |
| 3 Bladder   |                    |
| 4 External iliac vein                                   |                    |
| 5 External iliac artery                                 |                    |
| 6 Left testicular vein                                  |                    |
| 7 Left testicular artery                                |                    |
| 8 Obliquus externus abdominis muscle                    |                    |
| 9 Inguinal [abdominal] canal                            |                    |
| 10 Aponeurosis of obliquus externus abdominis muscle    |                    |
| 11 Cremaster muscle                                     |                    |
| 12 Internal spermatic fascia                            |                    |
| 13 External spermatic fascia                            |                    |
| 14 Cremasteric fascia                                   |                    |
| 15 Skin   |                    |
| 16 Scrotum  |                    |
| 17 Glans of penis                                       |                    |
| 18 Superficial fascia of penis                          |                    |
| 19 Deep fascia of penis                                 |                    |
| 20 Pubococcygeus muscle                                 |                    |
| 21 Deferent duct [spermatic duct]                       |                    |
| 22 Prostate   |                    |
| 23 Pubic bone   |                    |
| 24 Levator ani muscle                                   |                    |
| 25 Sphincter ani externus muscle                        |                    |
| 26 Anus [anal orifice]                                  |                    |
| 27 Seminal vesicle                                      |                    |
| 28 Ureter   |                    |
| 29 Rectal ampulla                                       |                    |
| 30 Internal urethral orifice                            |                    |
| 31 Medial umbilical fold                                |                    |
| 32 Fundiform ligament of penis                          |                    |
| 33 Cavernous body of penis                              |                    |
| 34 External opening of masculine urethra                |                    |
| 35 Navicular fossa of urethra [great lacuna of urethra] |                    |
| 36 Spongy part of masculine urethra                     |                    |
| 37 Perineal membrane                                    |                    |
| 38 Transversus perinei profundus muscle                 |                    |
| 39 Sphincter ani internus muscle                        |                    |
| 40 Rectovesical excavation                              |                    |
| 41 Prostatic utricle                                    |                    |
| 42 Anal canal   |                    |
| 43 Raphe of perineum                                    |                    |
| 44 Spongy body of penis                                 |                    |
| 45 Prepuce of penis                                     |                    |
| 46 Superficial dorsal vein of penis                     |                    |
| 47 White line of abdomen                                |                    |

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